

APPLICATION NOTE

Quantification of mitochondrial membrane potential

Using Celloger® Pro

■ Introduction

Mitochondria serve as the primary organelles responsible for generating energy within cells. In particular, Mitochondrial Membrane Potential (MMP) is considered a key indicator of mitochondrial function, as it is vital for ATP production.¹ Hence, MMP measurement has long been used as a marker for evaluating apoptosis, cell viability, and overall cellular health. Moreover, given that mitochondrial function is recognized to play a crucial role in conditions such as cancer, diabetes, and brain disorders like Alzheimer's and Parkinson's diseases², understanding and monitoring MMP become more significant.

Today, numerous fluorescent probes have been developed to visualize mitochondrial potential, including TMRE, TMRM, JC-1, JC-10, Mito-ID, MitoTracker and so on. Among these, JC-1 and JC-10 dyes exhibit changes in fluorescence wavelength depending on the probe's form. In healthy cells with normal MMP, JC-10 accumulates within the mitochondria, forming aggregates that emit red fluorescent wavelengths. On the other hand, in unhealthy cells with declined MMP, such as apoptotic cells, JC-10 scatters across the cytoplasm and converts into monomers, emitting green fluorescent wavelengths.³ In this application note, we introduce a procedure for evaluating mitochondrial membrane potential using Celloger® Pro.

■ Method

HeLa cells were seeded at a density of 3.5×10^4 cells per well in a 24-well plate. Following overnight incubation, the cells were treated with carbonyl cyanide 3-chlorophenylhydrazone (CCCP) (Sigma, C2759-100MG) at concentrations of 0.5, 1, 2, 4, 6, 8, 10 and 20 μM for 20 minutes. After a wash with HBSS buffer, the cells were stained with 10 μM JC-10 (AAT bioquest, 22204) solution in the dark for an hour. Finally, images were taken with Celloger® Pro, automated live cell imaging system with a 2X lens.

■ Result

We observed changes in fluorescence resulting from depolarization using JC-10 dye and carbonyl cyanide 3-chlorophenylhydrazone (CCCP), which acts as an uncoupler of oxidative phosphorylation. For quantitative comparison, we captured green and red fluorescence images using the Celloger® Pro with 2X lens for a wide field-of-view. And then, we employed analysis software to measure the intensity changes. As the concentration of CCCP increased, the membrane potential depolarized, leading to a decrease in red intensity and an increase in green intensity. (Fig. 1A). This trend persisted until a concentration of 20 μM was reached, at which point red intensity became almost negligible, causing the cells to appear entirely green. Figure 1B, captured with the Celloger® Pro using a 2X lens and cropped for magnification, enables the observation of red and green fluorescence distribution. The graph in Figure 1C, which depicts the green/red intensity ratio in relation to CCCP concentration as a dot plot, showed a gradual ascending pattern suitable for polynomial regression curve fitting.




Figure 1. Mitochondrial dysfunction depending on CCCP concentrations.

(A) The intensity of both green and red fluorescence varied with different CCCP concentrations. (B) In the control group, strong red dot-like signals were observed, whereas CCCP-treated cells displayed a prominent spread of green fluorescence throughout the cytoplasm. (C) The relationship between the green/red ratio and CCCP concentrations is modeled using a second-order polynomial curve ($R^2=0.97$).

■ Conclusion

Mitochondrial activity in cells is intricately regulated in response to intracellular conditions, including factors such as MMP, which can vary depending on the stage of the cell cycle.⁴ Therefore, it is important to analyze it within sufficiently wide areas to consider these variations, but collecting images from multiple regions is such a laborious task. Additionally, dyes like JC-10, which display two colors depending on their molecular forms, require individual detection using two different fluorescence channels. We successfully addressed these challenges by utilizing the Celloger® Pro, which offers features such as dual color fluorescence and bright-field imaging. Moreover, it efficiently captures images from multiple positions by automatically moving the camera.

Reference

1. Zorova, Ljubava D., et al. "Mitochondrial membrane potential." *Analytical biochemistry* 552 (2018): 50-59.
2. Bazhin, Arkadiy A., et al. "A bioluminescent probe for longitudinal monitoring of mitochondrial membrane potential." *Nature chemical biology* 16.12 (2020): 1385-1393.
3. Sivandzade, Farzane, Aditya Bhalerao, and Luca Cucullo. "Analysis of the mitochondrial membrane potential using the cationic JC-1 dye as a sensitive fluorescent probe." *Bio-protocol* 9.1 (2019): e3128-e3128.
4. Hirusaki, Kotoe, et al. "Temporal depolarization of mitochondria during M phase." *Scientific Reports* 7.1 (2017): 16044.