

Accessories Multi-beads Shocker

Sample Holders and Accessories

■ YK-MSH004(S)UY Multi- Sample Holder















for 8 x 2ml/3ml pulverizing tubes, 4 x 22ml/50ml pulverizing tubes

● Operational Accessories included with YK-MSH004(S)UY Sample Holder

Pulverizing tube adapters			Pulverizing media (metal rods)		
Product code	Description		Product code	Description	
YK-ADP03F2UY 1 pc.	50ml pulverizing tube adapter for 3ml and 2ml pulverizing tubes. Inserted into the 50ml tube hole. 3ml tubes fit directly inside, but 2ml tubes need be fitted with a 2ml adapter. One or two samples can be pulverized simultaneously.		YK-MC0218R(S)UY 10 pcs.	metal rod used with 2ml pulverizing tubes	
YK-ADP023(S)UY 8 pcs.	for holding 2ml tubes and is inserted into the 3ml tube holes		YK-MC0316(S)UY 10 pcs.	metal rod used with 3ml pulverizing tubes	
YK-ADP10UY 4 pcs.	for holding 10ml tubes and is inserted into the 50ml tube holes		YK-MC1020RUY 4 pcs.	metal rod used with 10ml pulverizing tubes	
YK-ADP22UY 4 pcs.	for holding 22ml tubes and is inserted into the 50ml tube holes		YK-MC1025RUY 4 pcs.	metal rod used with 10ml pulverizing tubes	
Transparent pulverizing tubes (used w/ liquid nitrogen)			YK-MC2225RUY 4 pcs.	metal rod (short) used with 22ml pulverizing tubes	
Product code	Description		YK-MC2235RUY 4 pcs.	metal rod (long) used with 22ml pulverizing tubes	
YK-ST0220PC(S)UY 100 pcs.	2ml transparent PC tubes used for freeze pulverizing		YK-MC5028RUY 4 pcs.	metal rod (short) used with 50ml pulverizing tubes	
YK-ST0320PCFUY 100 pcs.	3ml transparent PC tubes used for freeze pulverizing		YK-MC5038RUY 4 pcs.	metal rod used (long) with 50ml pulverizing tubes	
YK-ST1010PCUY 10 pcs.	10ml transparent PC tubes used for freeze pulverizing		YK-MC5055RUY 4 pcs.	metal rod used with 22ml and 50ml pulverizing tubes	
YK-ST2210PCUY 10 pcs.	22ml transparent PC tubes used for freeze pulverizing				
YK-ST5010PCRUY 10 pcs.	50ml transparent PC tubes used for freeze pulverizing				

Sample Holders and Accessories

● Included Accessories with YK-MSH004(S)UY Sample Holder

Opaque pulverizing tubes (used w/ liquid nitrogen)			Magnets		
Product code	Description		Product code	Description	
YK-ST0220FC(S)UY 100 pcs.	2ml opaque PE pulverizing tubes		YK-MG136UY 1 pc.	small (for removal of metal rod)	
YK-ST0320FCUY 100 pcs.	3ml opaque PE pulverizing tubes		YK-MG2810UY 1 pc.	large (for removal of metal rod)	
YK-ST1010CUY 10 pcs.	10ml opaque PE pulverizing tubes				
YK-ST2210CUY 10 pcs.	22ml opaque PE pulverizing tubes				
YK-ST5010CUY 10 pcs.	50ml opaque PE pulverizing tubes				
Translucent pulverizing tubes used within an approximate temperature range of -10°C ~ 20°C.					
Product code	Description				
YK-ST0250F(S)UY 100 pcs.	2ml translucent PP pulverizing tubes				
YK-ST0350FOUY 100 pcs.	3ml translucent PP pulverizing tubes				
YK-ST1010PPUY 10 pcs.	10ml translucent PP pulverizing tubes				
YK-ST2210PPUY 10 pcs.	22ml translucent PP pulverizing tubes				
YK-ST5010PPRUY 10 pcs.	50ml translucent PP pulverizing tubes				

Accessories Multi-beads Shocker


Sample Holders and Accessories

■ YK-SHAL224F(S)UY (2ml x 24) Sample Holder




● Included Accessories with YK-SHAL224F(S)UY Sample Holder


Pulverizing media (metal rod)

Product code	Description	
YK-MC0218R(S)UY 10 pcs.	metal rod used with 2ml pulverizing tubes	


Transparent pulverizing tube (used w/ liquid nitrogen)

Product code	Description	
YK-ST0220PC(S)UY 100 pcs.	2ml transparent PC tubes used for freeze pulverizing	



Opaque pulverizing tube (used with liquid nitrogen)

Product code	Description	
YK-ST0220FC(S)UY 100 pcs.	2ml opaque PE pulverizing tubes	

Translucent pulverizing tube used within an approximate temperature range of -10°C ~ 20°C.

Product code	Description	
YK-ST0250F(S)UY 100 pcs.	2ml translucent PP pulverizing tubes	

Magnets

Product code	Description	
YK-MG136UY 1 pc.	small (for removal of metal rod)	
YK-MG2810UY 1 pc.	large (for removal of metal rod)	


Sample Holders and Accessories

■ YK-SH1002(S)UY (100ml x 2) Sample Holder






● Included Accessories with YK-SH1002(S)UY Sample Holder




Pulverizing media (metal rod)

Product code	Description	
YK-MC100R(S)UY 2 pcs.	metal rod used with 100ml pulverizing tubes	


Transparent pulverizing tubes (used w/ liquid nitrogen)

Product code	Description	
YK-ST10010PCR(S)UY 10 pcs.	100ml PC tubes used for freeze pulverizing with inner caps and fitted silicon O-rings	 
YK-CP10010PCUY 10 pcs.	Outer cap for 100ml PC pulverizing tubes	

Translucent pulverizing tubes used within an approximate temperature range of -10°C ~ 20°C.

Product code	Description	
YK-ST10010PPR(S)UY 10 pcs.	100ml pulverizing PP tube with inner cap	 
YK-CP10010PPUY 10 pcs.	100ml pulverizing PP tube outer cap	

Magnets

Product code	Description	
YK-MG2810UY 1 pc.	large (for removal of metal rod)	

Accessories Multi-beads Shocker


Sample Holders and Accessories

■ YK-SHAL2318(S)UY (2 - 3ml x 18) Sample Holder





● Included Accessories for YK-SHAL2318(S)UY Sample Holder



Pulverizing tube adapter

Product code	Description	
YK-ADP023(S)UY 18 pcs.	For holding 2ml tubes and is inserted into the 3ml sample holder tube holes	



Transparent pulverizing tubes (used w/ liquid nitrogen)

Product code	Description	
YK-ST0220PC(S)UY 100 pcs.	2ml transparent PC tubes used for freeze pulverizing	
YK-ST0320PCFUY 100 pcs.	3ml transparent PC tubes used for freeze pulverizing	



Opaque pulverizing tubes (used w/ liquid nitrogen)

Product code	Description	
YK-ST0220FC(S)UY 100 pcs.	2ml opaque PE pulverizing tubes	
YK-ST0320FCUY 100 pcs.	3ml opaque PE pulverizing tubes	



Pulverizing tubes used within an approximate temperature range of -10°C ~ 20°C.

Product code	Description	
YK-STC0220FSEUY 100 pcs.	2ml pulverizing PP tube with o-ring * can be used with metal rod as well as glass or zirconia beads	
YK-STC0350FOUY 100 pcs.	3ml pulverizing PP tube with o-ring * can be used with metal rod as well as glass or zirconia beads	

Pulverizing media (metal rod)

Product code	Description	
YK-MC0218R(S)UY 18 pcs.	metal rod used with 2ml pulverizing tubes	
YK-MC0316(S)UY 18 pcs.	metal rod used with 3ml pulverizing tubes	

Magnets

Product code	Description	
YK-MG136UY 1 pc.	small (for removal of metal rod)	
YK-MG2810UY 1 pc.	large (for removal of metal rod)	


Sample Holders and Accessories

■ YK-SHC224F(S)UY (2ml x 24) Sample Holder



● Included Accessories with YK-SHC224F(S)UY Sample Holder



Pulverizing media (metal rod)

Product code	Description	
YK-MC0218R(S)UY 10 pcs.	metal rod used with 2ml pulverizing tubes	



Pulverizing media (glass or zirconia beads)

Product code	Description	
YK-YGB01UY 500g	0.1mm glass beads	
YK-YGB05UY 500g	0.5mm glass beads	
YK-YZB01UY 200g	0.1mm zirconia beads	
YK-YZB05UY 200g	0.5mm zirconia beads	

Pulverizing tubes used within an approximate temperature range of -10°C ~ 20°C.

Product code	Description	
YK-STC0250FSUY 100 pcs.	2ml pulverizing PP tube used with pictured sample holder * use with glass or zirconia beads	
YK-STC0220FSE 100 pcs.	2ml pulverizing PP tube with o-ring * use with metal rod and glass or zirconia beads	

Magnets

Product code	Description	
YK-MG136UY 1 pc.	small (for removal of metal rod)	
YK-MG2810UY 1 pc.	large (for removal of metal rod)	

Accessories Multi-beads Shocker


Sample Holders and Accessories

■ YK-SHC2318(S)UY (2-3ml x 18) Sample Holder





● Included Accessories for YK-SHC2318(S)UY Sample Holder



Pulverizing tube adapter

Product code	Description	
YK-ADP023(S)UY 18 pcs.	For holding 2ml tubes and inserted into the 3ml sample holder tube holes	



Pulverizing media (metal rod)

Product code	Description	
YK-MC0218R(S)UY 18 pcs.	metal rod used with 2ml pulverizing tubes	
YK-MC0316(S)UY 18 pcs.	metal rod used with 3ml pulverizing tubes	


Pulverizing media (glass or zirconia beads)

Product code	Description	
YK-YGB01UY 500 mg	0.1mm glass beads	
YK-YGB05UY 500 mg	0.5mm glass beads	
YK-YZB01UY 200 mg	0.1mm zirconia beads	
YK-YZB05UY 200 mg	0.5mm zirconia beads	

Pulverizing tubes used within an approximate temperature range of -10°C ~ 20°C.

Product code	Description	
YK-STC0220FSEUY 100 pcs.	2ml pulverizing PP tube with o-ring * can be used with metal rod as well as glass or zirconia beads	
YK-STC0350FOUY 100 pcs.	3ml pulverizing PP tube with o-ring * can be used with metal cone as well as glass or zirconia beads	

Magnets

Product code	Description	
YK-MG136UY 1 pc.	small (for removal of metal rod)	
YK-MG2810UY 1 pc.	large (for removal of metal rod)	